

Scientific

Distributed by: ABQ Industrial LP USA Tel: +1 (281) 516-9292 / (888) 275-5772 eFax: +1 (866) 234-0451 Web: https://www.abqindustrial.net E-mail: info@abqindustrial.net

. 000

TAR STAR

DATA

LAQUA

Touchscreen precision. The new benchmark.

рН	ORP	ION	Conductivity			
Resistivity	Total Disso	Total Dissolved Solids				

NEAS

Benchtop Water Quality Meters

F-70/DS-70 Series



Ideas for you

Our concept originated from you

"Meters and electrodes become dirty so often, I wish I could keep them clean all the time."

"It would be great if I could quickly visualize the calibration and measurement results, as well as the status of the electrodes."

"If a problem occurs, I want it solved immediately!"

- "It' s a pain to have to look through the instruction manual."
- "Electrode stands are actually not that user-friendly."
- "I want the electrode stand to move freely according to the location and what I'm using it for."

"Robust with high precision."

"Which electrode actually matches my need?"

"I want stable measurements every day."

We listened closely to each of our customers' comments, and applied what we heard to our next generation analyzer. HORIBA is proud to announce LAQUA, the water quality analyzer that answers all of your needs.



LAQUA

Stress free operation / Smart navigation	• P03~05
Essence of technology New pH electrode	• P06~08
Product selection guide and packages meter + electrode	• P09~10

LAQUA Our promise | We want to respond flexibly to your analysis needs. Containing a word "water" in Latin, our commitment to provide everything you expect from a water quality analyzer is distilled in our new brand, LAQUA.

Intuitive and easy to use touch panel operation

Intuitive control with the large capacitive touch panel. Smart navigation provides step-by-step guidance for trouble-free operation. Easy to clean glass top and round body, LAQUA is both easy and fun to use!

SMART

Operation buttons are reduced to the bare minimum

CAL	MEAS	
Ъ	\sim	

Calibration

DATA

Measurement Data Management

Simply slide your finger across the screen to switch displays

Switch between digital, graphic, and analog displays during measurement with just the flick of a finger. No need for complex actions.

2-channel simultaneous measurement and display

pH value and a second measurement (such as ORP, ion, electrical conductivity) can be displayed simultaneously.

aQL

AS MEAS	A	÷ 2	2011/06/20 18:00
CH1 ►			25.0°C
pH 🕨		7.	000
CH2 ►		-	25.0°C
COND		0.	000 ms/m
STARTキーで測	『定を『	報告しま (ਰਾ
801		5	START





Accurate calibration for measurement precision.

Correct calibration is done under "stable" conditions. Calibration performed under unstable conditions is one of the big causes of measurement error. Calibration response is visualized as numerical data or a graph. With

LAQUA, you are sure about your calibration validity

Calibration Assistance Function

You can tell measurement value has stabilized when the graph has stabilized and the calibration stability values become smaller. "Stability" checking at a glance!



www.horiba.com/laqua

NAVIGATION

Functions for Validation

and Usability

For compatible models please

see the P14 body specification.

Enjoy hassle-free operation with on-screen settings confirmation, maintenance information, and troubleshooting tips

the WEB



 Customizable auto hold function for calibration and measurement
 Periodic inspection mode: JIS/Pharmacopeias/Digital Simulator
 Digital memory: Maximum 2,000 sets of measurement data can be recorded (999 sets for F-71/F-74BW/DS-71 models)
 Simultaneous connection to a GLP/GMP compatible printer and PC
 Customizable print function
 Save data onto a USB flash drive
 USB PC Communication: Data storage software available as a free download for registered users.
 Multi-language support (Japanese, English, Chinese, Korean)
 FDA21CFR Part 11 (Please ask for quotation)

04 | LAQUA

LAQUA's free arm electrode stand can handle any container size or position.

The stand-alone free arm electrode stand can be moved wherever you like, vertically or horizontally.

You can also use the long electrode stand* with a telescopic shaft when working with large beakers.

360°



With the long electrode stand*, you can prepare small quantities of standard solution for calibration or large capacity containers of buffer solution without having to detach and reattach the electrodes.

HORIBA electrode technology gives you the fusion of high accuracy and ease of use

ELECTRODE

HORIBA electrode is now even tougher and responds faster.



Enhanced stability and minimized drift

Integrating two new technologies for faster response times and optimal performance.



pH fast response glass membrane (U.S. Patent No. 8262877) The membrane contains HORIBA's unique combination of rare earth metals to improve response time by twofold and to increase

durability against chemical attack.





Reference electrode with increased stability (Patent pending)

Covering the internal electrode with a cation-conductive hollow fiber membrane, liquid junction clogging by silver ions and silver complex ions is reduced to 1/1000 of the conventional technology. Furthermore, maintained internal solution concentration ensures a stable standard electrical potential.

ToupH electrodes are now even stronger

HORIBA's glass membrane molding technology achieves strengths more than 10 times the Japanese Industrial Standards (strength tests).



New dome-shaped construction boosts strength in all directions!



New • pH electrode

AQUA

Not just "unbreakable." New flat sensor innovations allow the measurement of trace sample droplets or the measurement of solid sample surface.



What is an ISFET(semiconductor sensor)?

ISFET is the abbreviation of Ion Sensitive Field Effect Transistor. The response membrane is equipped with semiconductor based sensor.

ISFET features

1.Will not crack or break like conventional glass electrodes

2.The sensor is flat and very small in size, enabling the measurement of extremely small samples
3.Easy handling and maintenance – simply clean with a toothbrush
4.Can be stored dry

The flat electrode has less than 100 $\!\mu\text{m}$ distance between the housing and the sensor.

The unique structure enables to measure miniscule amount of moisture on the surface of solid objects and prevents bubbles from trapping on the sensor when measuring samples in a beaker.

Effects of static electricity is reduced

The combination of HORIBA's unique semiconductor device construction and improved static protection circuit means that the effects of static electricity, once the Achilles heel of semiconductor sensors, are greatly reduced.

Precision—pH electrodes from HORIBA which answers your needs.

Stable measurement for a wide range of samples. Standard **ToupH** electrode (9615-10D)

STANDARD ToupH 📾 📾

High stability and drift reduction. No more worries about the timing of your measurement value readings.

- Uses responsive glass that is 10 times stronger than JIS standard. The domed shape provides strength in all directions, greatly reducing damage concerns.
- Constructed with smooth surfaces for easy wiping and cleaning.

Recommended

Perfect for preparing buffers. Can be used on a wide range of aqueous test solutions.

For extremely small samples Micro (ToupH) electrode (9618-10D)





This pH electrode with temperature compensation sensor can take measurements from samples as small as 50µL, the smallest in the world.

- Our original manufacturing technology (Japanese Patent No. 4054245) is used to produce 2-ply piping 3mm in diameter.
- Compatible with extremely small containers such as micro tubes etc.
- The temperature sensor is located at the tip for high-speed temperature response. Refrigerated samples can be measured without needing to wait for them to return to room temperature.



LAQUA

LAQUA

IAQUA

Recommended

Can be used for a wide range of aqueous solutions, including those that cannot be obtained in large quantities. We recommend using our specialized cleaning solution after measuring samples that contain proteins.

For using a large container Long **ToupH** electrode (9680-10D)



283 mm length & 8 mm diameter. The long, thin design makes this electrode perfect for measuring in large containers and test tubes.

 Uses responsive glass that is 10 times stronger than JIS standard. The domed shape provides strength in all directions, greatly reducing damage concerns.

Recommended

For measuring samples such as microbe culture fluids in test tubes. We recommend that it be used with the long type electrode stand (FA-70L).

For highly viscous samples Sleeve **ToupH** electrode (9681-10D)



SLEEVE Touph





Stable measurement can also be achieved for high viscous samples.

• The liquid junction section is constructed with a moveable sleeve that can be rinsed clean, preventing highly viscous samples from clogging the liquid junction, and maintaining stable measurement performance

Recommended \

For highly viscous samples and solutions, and samples that contain non-aqueous solvents (such as cosmetics or paints). We recommend that you take measurements while using the graph display function to confirm stable responses. (We recommend washing with a neutral detergent after use with samples that contain oil.)

For the surface of solid samples Flat ISFET pH electrode (0040-10D)







LAQUA

The sensor is located on the flat surface of the electrode tip, with less than a 100 μm protrusion from the housing.

- Measurements can be made from a minute amount of moisture on the solid sample surface.
- Use of a semiconductor sensor means there are no concerns that the electrode will be damaged.
 - Also perfect for measuring samples in shallow containers such as Petri dishes.
 - Repalceable sensor

Recommended \

For surface measurement of gelatinous materials such as nutrient agar, and food samples such as meat. Evaluation of sheet materials such as cloth or paper. (If the sample only has a small amount of moisture, pure water etc. is required. We recommend washing with a neutral detergent after use with samples that contain oil.)

For easy and safe measurement inside solid samples (0030-10D)

NEEDLE ISFET



The sharp tip can pierce solid

material to take measurement within the sample.

- Use of a semiconductor sensor means there are no concerns that the electrode will be damaged.
- Repalceable sensor

Recommended \

For measuring inside foodstuffs, such as fruits, vegetables and bread. (We recommend washing with a neutral detergent after use with samples that contain oil.)

For stable measurement of tap water Low conductivity/Low buffer capacity pH electrode (9630-10D)

For TAP WATER 📾 📾



Using the high-purity glass membrane, faster stable measurement is possible at a low electrical conductivity and low buffer capacity sample

 It enables the measurement within 90 seconds measurement (Auto hold) for tap water by using the conditioning liquid(model name 230). (95% Response is within 60 seconds)

Recommended \

It is ideal for water quality testing in the water purification plant.

LAQUA



Recommended Packages

Complete sets with meter, electrode, and standard solutions



Custom LCD





Benchtop pH / ORP Custom LCD Meter Set, complete with

- electrode stand
- glass pH electrode (9615-10D)
- pH 4/7/10 buffers & 3.33M KCl referece electrolyte (502-S)

DS-71A-S

Benchtop Conductivity / Resistivity / Salinity / TDS Custom LCD meter, complete with

- electrode stand
- conductivity electrode (3552-10D)
- conductivity standard solutions 84 uS, 1413 uS, 12.88 mS & 111.9 mS (503-S)



F-74BW-A-S

Benchtop pH / ORP / ION / Conductivity / Resistivity / Salinity / TDS Custom LCD meter complete with

- electrode stand
- glass pH electrode (9615-10D)
- pH 4/7/10 buffers & 333M KCl referece electrolyte (502-S)
- conductivity electrode (3552-10D)
- conductivity standard solutions 84 uS, 1413 uS, 12.88 mS & 111.9 mS (503-S)

Touch Screen Color LCD

2000



F-72A-S

Benchtop pH / ORP / ION Color Touch Screen Meter Set, complete with

- electrode stand
- glass pH electrode (9615-10D)
- pH 4/7/10 buffers & 3.33M KCl referece electrolyte (502-S)

F-73A-S

Benchtop pH / ORP / ION, Dual Channel Color Touch Screen Meter set,

- complete with
- electrode stand
- glass pH electrode (9615-10D)
- pH 4/7/10 buffers & 3.33M KCl referece electrolyte (502-S)



Benchtop Conductivity / Resistivity / Salinity / TDS Color Touch Screen meter, complete with

- electrode stand
- conductivity electrode (3552-10D)
- conductivity standard solutions 84 uS, 1413 uS, 12.88 mS & 111.9 mS (503-S)

F-74A-S

Benchtop pH / ORP / ION / Conductivity / Resistivity / Salinity / TDS Color Touch Screen meter complete with

- electrode stand
- glass pH electrode (9615-10D)
- pH 4/7/10 buffers & 3.33M KCI referece electrolyte (502-S)
- conductivity electrode (3552-10D)
- conductivity standard solutions 84 uS, 1413 uS, 12.88 mS & 111.9 mS (503-S)

pH Electrode Selection Guide

				3-in-1 ELECTRODES (ToupH)				ISFET ELE	ISFET ELECTRODES		3-in-1 ELECTRODES	
			PLASTIC	STANDARD ToupH	LONG ToupH	MICR0 ToupH	SLEEVE ToupH	NEEDLE ISFET	FLAT ISFET	SLEEVE	NON- AQUEOUS	
		9625-10D	9615-10D	9680-10D	9618-10D	9681-10D	0030-10D	0040-10D	6367-10D	6377-10D		
Applicable temperature range (°C)		0-100	0-100	0-100	0-60	0-60	0-60	0-60	0-60	0-60		
Specification	Diameter (mm)		16	12	8	3	12	15	10	12	12	
opeemeation	Position of liquid junction (approx. mm)		15	13	21	6	26	11	0.1	10	23	
	Length (mm)		150	198	283	185	203	190	190	150	150	
pH - Sample	e Conditions											
		Normal (over 100 mS/m)	۲	۲	۲	۲	۲	۲	۲	۲	۲	
	Conductivity	Low (approx.10 ~100 mS/m					0				۲	

	Conductivity	Low (approx.10 ~100 mS/m					0				•
		Very low (approx. 5 ~100 mS/m					0				۲
		High (approx. 5 S/m)	0	0	0		۲				
Aqueous Solution	Strong alkaline (pH 10-12)			0	0		0			0	
	Strong acidity (pH 0-2) * Except HF sample			۲							
	Quick heat change (within 50°C)		۲								
	High viscosity (approx. 5 Pa·S)						۲			0	۲
	Containing non-aqueous solvent			0	0	0	0	0	0	0	۲
	Suspension			0	0	0	۲	0	0		۲
Solid/Semisolid	Inside							۲			
SUIIU/Semisuiiu	Surface								۲		

pH - Sample Conditions

pri - Sample											
	Microtube/plate	(> 50 µL)	×	×	×	۲	×	×	×	×	×
	NMR tube	ø5 mm ID > ø4 mm	×	×	×	×	×	×	×	×	×
	Ampule	>ø4 mm				۲					
Comple	Micro container (> 2 mL)				0	۲					
Containers	Tube	ID:13 mm, L:100 ~ 150 mm			۲						
Containers	Beaker	10 mL ~ 1 L	۲	۲	0	0	0	0	0	0	0
	Large container (> 1 L)		0	0	۲						
	Petri dish								۲		
	Droplet		×	×	×	×	×	×	۲	×	×

pH - Typica	al Samples									
	Pure/ion-exchange water (approx. 0.1 mS/m)									۲
	Distilled water (approx. 0.5 mS/m)		0							۲
Water	Tap/drinking water (approx. 10 mS/m)	0	0			0				۲
water	Surface water		0			0				۲
	Pharmaceutical water		0			0				0
	Enviromental water/acid rain	0	0			0				0
	Caustic/strong acid (Except HF sample)		۲			0				
	Hydrofluoric acid									
Chemical	Organic solvent	×					×	×		0
reagent/solvent	KCI-reactive solution	×	×	×	×	×	×	×	×	×
rougont, contoint	Surfactant		0			۲				0
	Water-based paint		0			۲				0
	Dye/coloring agent					۲				0
	Protein-containing sample		0		0	۲			0	
	Medicinal preparation				0	0				0
Pharmaceutical/ biology sample	Enzyme solution			0	۲					
	Tris buffer		۲		0	0				
	Suspension		0			۲				۲
	Agar medium							۲		
	Jam		0			۲	(inside)	(surface)		0
	Meat/fish						(inside)	(surface)		
	Fruit/vegetable						(inside)	 (surface) 		
Food	Dough						(inside)	(surface)		
	Honey						(inside)	(surface)		۲
	Cheese/butter						(inside)	(surface)		
	Yogurt	0	0			0	(inside)	 (surface) 	0	
	Beer	0	0			۲			0	۲
Beverage/	Milk		0			۲			0	0
seasoning	Carbonated drink/juice/sauce/soy sauce		0			۲			0	0
	Mayonnaise/ketchup		0			۲				0
	Beauty cream/mascara		0			۲	0			0
Cosmetic/	Gel/soap/shampoo		0			۲				0
lotion	Hairdye lotion		0			۲				0
	Emulsified liquid		0			0				۲

O Recommended \bigcirc Can be measured $~\times$ Prohibited or risk of damage

Electrodes/Accessories



pH Electrode					
	Description	Model	Temp. range (°C)	pH range	Part No.
	Plastic body	9625-10D	0~100*1	0~14	3200360505
	Standard ToupH	9615-10D	0~100	0~14	3200366539
	Sleeve ToupH	9681-10D	0~60	0~14	3200366572
	Long ToupH	9680-10D	0~100*1	0~14	3200366560
	Micro ToupH	9618-10D	0~60	0~14	3200366552
Combination (3-in-1)	Sleeve	6367-10D	0~60	0~14	3014079136
pH electrode	For measurement of low-conductivity water and non-aqueous solvents	6377-10D	0~60	0~14	3014093085
	Needle type	6252-10D	0~60	0~12	3014080850
	For Tap water	9630-10D	0~100	0~14	3200528726
	For Hydrofluoric acid sample	9631-10D	0~60	2~12	3200524119
	For Strong alkali sample	9632-10D	0~100	0~14	3200524120
	Needle type ISFET	0030-10D	0~60	0~14	3014028323
ISFET	Flat type ISFET	0040-10D	0~60	0~14	3200367925
pH electrode	Needle type ISFET(0030-10D) sensor	0131	0~60	0~14	3014028400
	Flat type ISFET(0040-10D) sensor	0141	0~60	0~14	3200367926
Combination	For very slender test tubes	6069-10C	0~60	0~14	3014081107
pH electrode	Flat type	6261-10C	0~50	0~12	3014081807
	Standard type	1066A-10C	0~100	0~14	3014080432
Glass pH electrode	For measurement of low-conductivity water and non-aqueous solvents.	1076A-10C	0~100	0~14	3014093084
Deference electrode	Standard type	2060A-10T	0~100	—	3014080434
Reference electrode	Double-junction type	2565A-10T	0~100	—	3014080436
Temperature electrode	For temperature compensation and measurement	4163-10T	0~100	_	3014080375
ORP electrode	Platinum 3-in-1 type	9300-10D	0~ 60		3014046710

Conduc	ctivity Cell		 Conductive material: Platinum ring: Body housing: Glass except 9382- 	s coated with platinum IOD - Plastic	black		
Cell constant cm ⁻¹ (m ⁻¹)		Model	Range cm-1(m-1)	Minimum Volume (mL)	Application	Temp. range (°C)	Part No.
	0.1 (10)	3551-10D	0.1 µS~10 mS (10 µS~1 S)	50	For low conductivity water (deionized water or other)	0~60	3014081712
Immersion type	1 (100)	9382-10D	1 µS~100 mS (0.1 mS~10 S)	20~30	Waterproof. For general purposes	0~80	3014046709
	1 (100)	3552-10D	1 µS~100 mS (0.1 mS~10 S)	15	For general purposes	0~100	3014081545
	10 (1000)	3553-10D	10 µS~1 S (1 mS~100 S)	50	For high conductivity water	0~60	3014081714
	0.1 (10)	3561-10D	0.1 μS~10 mS (10 μS~1 S)	10	For low conductivity water (pure water or other)	0~60	3014082350
Flow turns	1 (100)	3562-10D	1 µS~100 mS (0.1 mS~10 S)	16	For general purposes	0~60	3014082513
Flow type	10 (1000)	3573-10C	10 µS~1 S (1 mS~100 S)	4	For high conductivity water	0~60	3014082590
	10 (1000)	3574-10C	10 µS~100 mS (1 mS~10 S)	0.25	For column chromatography using a very small amount of sample	0~60	3014082592

Ion Selective Electrode		All ion electrodes (except combination el Please be aware of the hindering ion and	ctrodes only.	Replacement Tip		
Electrode name	Model	Measuring range	Interfering ion influence"	Part No.	Model	Part No.
Combination Chloride ion electrode*	6560-10C	0.4~35,000 mg/L Cl⁻	Br ⁻ =0.03 NO ₃ , F ⁻ , HCO ₃ ⁻ , SO ₄ ²⁻ , PO ₄ ²⁻ =1,000	3014093430	7660	3014093436
Combination Fluoride ion electrode*	6561-10C	0.02~19,000 mg/L F	(ex. Al ³⁺ , Fe ³⁺)coexisted and foamed the complex.	3014093431	7661	3014093438
Combination Nitrate ion electrode*	6581-10C	0.62~62,000 mg/L NO₃⁻	CH3COO ⁻ =300 SO4 ²⁻ =Over 1000	3014093432	7681	3014068364
Combination Potassium ion electrode*	6582-10C	0.04~39,000 mg/L K*	Li ⁺ , Na ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ =Over 1000	3014093433	7682	3014069795
Combination Calcium ion electrode*	6583-10C	0.4~40,080 mg/L Ca ²⁺	Mn ²⁺ =500 Mg ²⁺ =1,000 Na ⁺ , K ⁺ , Ba ²⁺ , NH ₄ ⁺ =Over 1,000	3014093434	7683	3014068795
Combination Ammonia electrode*	5002A-10C	0.1~1,000 mg/L NH₃	_	3014093560	membrane (NH₃)	3014067083

*1 The selection coefficient is a ratio of the limit concentration of coexisting ions (mol/L) to the ion concentration to be measured (mol/L); A value of 1000 means that the coexisting ions can be permitted up to 1000 times the ion measured and "N/A" means that chemical change occurs in the solid response membrane.

pH Solution Kits										
Name	Туре	Specification	Volume	Part No.						
NIST pH Buffer Solution Kit	501-S	(4.01/6.86/9.18/KCI Reference)	250mL ea	3999960015						
USA pH Buffer Solution Kit	502-S	(4.01/7.00/10.01/KCI Reference)	250mL ea	3999960016						
pH Solutions										
	500-2	pH 1.68	500ml	3999960028						
	500-4	pH 4.01	500ml	3999960029						
	500-686	pH 6.86	500ml	3999960030						
Buffer Solution at 25°C	500-7	рН 7.00	500ml	3999960031						
	500-9	рН 9.18	500ml	3999960032						
	500-10	pH 10.01	500ml	3999960033						
	500-12	pH 12.46	500ml	3999960034						

Conductivity Solution Kit							
Name	Туре	Specification	Volume	Part No.			
Conductivity Standard Solution Kit	503-S	(84 uS/1413 uS/12.88 mS/111.8 mS)	250ml ea	3999960017			
Conductivity Solutions							
Conductivity Standard Solution at 25°C	500-21	84 uS	500ml	3999960035			
	500-22	1413 uS	500ml	3999960036			
	500-23	12.88 mS	500ml	3999960037			
	500-24	111.8 mS	500ml	3999960038			

ORP				
Name	Туре	Specification	Part No.	
Powder for ORP Standard Solution	160-51	89 mV For 250 mL (10 packets per set)	3200043618	
	160-22	258 mV For 250 mL (10 packets per set)	3200043617	

Internal Filling Solution for Electrodes					
Name	Туре	Specification	Volume	Part No.	
Internal Filling Solution for pH Combination Electrode	525-3	3.33 M KCI	250ml	3999960023	
Internal Filling Solution for Reference Electrode	300	3.33 M KCI	250ml	3200043640	

Accessories					
		Name	Part No.		
		Printer (for GLP/GMP compliance) Cable sold separately, Plain paper	3014030147 (230v) 3014030146 (120v)		
Printer	Printer Printer cable	Printer cable (1.5 m)	3014030148		
		Printer paper (20 rolls)	3014030149		
	Ink ribbon Printer paper	Ink ribbon (5 pcs/set)	3014030150		
Power	AC adapter	AC adapter cable set for LAQUA	3014031952 (230v)		
	, and the datapier	meters. (AC adaptor 1.6 m, cable 1 m)	3014031951 (120v)		
For Inspection	Tan.	Digital simulator X-51 (pH, mV, ION, DO simulator)	3014028368		
For inspection	X-51 X-52	Digital simulator X-52 (Conductivity simulator)	3014028370		
Motor	\cap	LCD protection sheet (2 pcs/pack)	3200382462		
Accessories	LCD protection sheet cover	Protection cover (Protects the meter for F-70, DS-70 series)	3200382441		
		USB cable (Cable to connect meter and PC.)	3200373941		
Communication and Output USB cable Serial	\bigcirc	Analog cable (Analog (alarm) output cable)	3014030152		
	USB cable Serial cable	Serial cable (Cable to connect meter and PC (Serial, 9 pins))	3014030151		
Electrode Stand (images on the right)		FA-70S Electrode stand (adjustable type) (Free-standing type. Height 384 mm)	3200382557		
		FA-70L Electrode stand (long type) (Free-standing type. Height 450~650mm)	3200382560		
	Arm for electrode stand	Arm for electrode stand (For FA-70S, FA-70L)	3200373991		
Electrode Accessories		Sensor Holder (Used for Mounting Electrode Stand, 2 pcs.)	3200373961		
		Electrode Protection Cap (Standard) (For 9615-10, 9618-10D, 9681-10D pH Electrode, 3 pcs.)	3200382477		
		Electrode Protection Cap (Standard) (For 9621-10D, 9625-10D, 9630-10D, 9631-10D, 9632-10D, 6367-10D, 6377-10D, 6252-10D, 6261-10C, 10664-10C, 1076-10C, 2060-10T, 9300-10D, 9382-10D, 3552-10D pH Electrode, 5 pcs.)	3200043508		
		Electrode Protection Cap for Long Electrode (For 9678/9680 pH Electrode, 1 pc.)	3200382482		





Standard Electrode Stand FA-70S (384mm)

When arm is housed 450mm

650mm

Long Type Electrode Stand FA-70L (450~650mm)

-

LAQUA F-70/DS-70 series specifications

		F-71	F-72	F-73	F-74	F-74BW	DS-71	DS-72	
	Measurement method		Gla	ass electorode me	thod		—	—	
	Measurement range	рН 0.000~14.000 —						-	
	Display range	pH -2.000~19.999		pH -2.000~20.00	0	pH -2.000~19.999			
	Resolution	0.001 pH		0.01/0.001 pH	-	0.001 pH	-	_	
рН	Auto range select	-	•	•	•	_	_		
	Repeatability	±0.005 pH±1 digit		±0.001 pH±1 digi	t	±0.005 pH±1 digit			
	Repeatability check	5		5		5			
	Alarm limit of calibration	•	•	•	•	•	_	_	
	Periodical check	_	•	•	•	-	_	-	
	Measurement range	I		±1999.9 mV			_	-	
mV (ORP)	Resolution			0.1 mV			—	—	
	Repeatability	±0.1 mV±1 digit — —							
	Measurement range	0.0~100.0°C (-30.0~130.0°C)							
Iemperature	Resolution	0.1°C							
	Measurement method			lon cloatro					
	Measurement range	_		0.00 ug/l ~9	99 a/L (mol/L)				
	Resolution	_		3 signific	nificant digits			-	
ION	Repeatability	_		±0.5%F.	S.±1 digit		-	-	
1011	Periodical check	_	•	•	•	-	_		
	Calibration curve point	_	5	5	5	5	—	—	
	Addition method measurement	_	•	•	•	-	_		
	Measurement method	_	_	-		2 AC bipola	r method		
	Measurement range (Display range)	_	—	-	0.0 µS/cm~19.99 µS/cm : Cell constant 0.1/cm 0.000 mS/cm~199.9 mS/cm : Cell constant 1.0/cm 0.00 mS/cm~1999.0 mS/cm : Cell constant 10.0/cm				
	Resolution	_	_	_		0.05% of fi	ull scale		
Conductivity	Repeatability	_	_	_		±0.5%F.S.	±1 digit		
	Measurement unit selection	_	—	-	•	•	•	•	
	Distilled water temperature conversion	_	_	_	•	•	•	•	
	IP/EP/LISP/CP Pharmaceutical water anlication	_		_		_	_		
	Measurement method	_	_	_	•	Conversion from c	onductivity value		
	Measurement range (Display range)	_	_	-		0.00~80.00 ppt (0	.000%~8.000%))	
Salinity	Resolution	_	_	-		0.01 ppt (0	.001%)	,	
	Salt concentration calibration	_	_	-	•	•	•	•	
	Measurement method	_	_	-		Conversion from c	onductivity value	\$	
Resistivity	Measurement range (Display range)	_	_	-	0.0 Ω • cm~199.9 MΩ • cm : Cell constant - 0.1/cm 0.00 Ω • cm~19.99 MΩ • cm : Cell constant - 1.0/cm				
	Resolution	—	_	-		0.05%	F.S.		
	Repeatability	—	_	-	±0.5%F.S.±1 digit				
	Measurement method	_	-	-	Conversion from conductivity value (EN27888 or TDS Factor)			r TDS Factor)	
TDS	Measurement range (Display range)	_	_	-	0.01 mg/L~1000 g/L	. 0.01 mg/L^	~100 g/L	0.01 mg/L~1000 g/L	
	Resolution	- 1	1	2	2	0.01 m	1g/L	1	
	USB peripherals (Communication with PC) ¹	1	•	2	2	2			
Input/	USB host (USB memory)	-	•	•	•	-	_		
output	RS-232C (Printer/PC)	•	•	•	•	•	•	•	
	Analog output	_	٠	•	•	_	_	•	
	Memory number	999	2000	2000	2000	999	999	2000	
Data	Interval memory	•	•	•	•	•	•	•	
	ID input	•	•	•	•	•	•	•	
	Data search	_	•	•	•	-	_		
	Display	Custom LCD	Color graphic	LCD with capaciti	ve Touch Panel	Custom	LCD	color graphic LCD with capacitive Touch Panel	
Display	Dual component display	_	_	•	•	•	_	_	
1 5		_			- //(Japanese/English/	
	Multilanguage display		Japane	se/English/Chines	e/Korean	_		Chinese/Korean	
	Navigation function	_	•	•	•	-	-	•	
	User guide	-	•	•	•	_	-	•	
	Graph display	_	•	•	•	-		•	
	Printer connectivity (GLP/GMP)	•	•	•	•	•	•		
							-		
Function	AutoHold function	•	•	•	•	•	•	•	
	AutoHold setting	_	٠	•	•	-	-	•	
	Stability function (pH/ION)	-	•	•	•	-	-	•	
	Operator ID	—	٠	•	•	-	-	•	
	Security (password)	•	•	•	•	•	•	•	
	Version up function	•	•	•	•	•	•	•	
Ambient ter	nperature								
I Uwoi AC adaptor I UV ~ 240 V SU/00 HZ Dimensions 170 /W/\v174 /D\v73 /H\mm (Evaluation algobrade stand and AC adaptor)									
Power cons	sumption	Approx. 0.7 VA	170 (11)	AV 8 9. xorgq		Annroy	0.7 VA	Approx, 9.8 VA	
Mass of main unit		Approx. 500 g Approx. 700 g			Approx. 500 g Approx. 700 g				

*1 USB cable sold separately. Software can be download by web registration.



AQUA

LAQUA .

Visit HORIBA's website!

Water Quality Analyzers www.horiba-water.com

With over 60 years of engineering excellence, HORIBA's diverse range of water quality analyzers and electrodes are ideal for everyday laboratory needs through to the most demanding of applications. Visit our website for a wealth of useful information and water quality measurement tips to help you obtain the best results in your work.



Benchtop Meters

Developed using extensive feedback from users, our new LAQUA meters deliver the best solution for water quality analysis. Our LAQUA website features an online 'Selection Guide' to enable you to find the perfect LAQUA meter and electrode for your need.

Handheld Meters

In the lab, in the field or anywhere you need it. LAQUA Handheld meters are designed for use with one hand and with an IP67 waterproof rating and shock-resistant casing. Meters can be used for long periods, even in dark places, making it ideal for field measurements in rivers and lakes.

Pocket Meters

Analyzing water quality is simplified when using our LAQUAtwin range of meters. Designed to produce accurate and reliable results. Anyone, anywhere, at any time can measure samples easily with a LAQUAtwin meter. See just how good they are at our website.

Electrodes

Various electrodes to match any application. A wide range of products for both benchtop and portable systems are available, including easy and reliable standard models, application-focused models for small samples or large containers, and special electrodes for specific sample characteristics.

SUPPORT HORIBA CUSTOMER SUPPORT SYSTEM

HORIBA offers a variety of services to conform to quality standards and international guidelines such as GLP, GMP and ISO

Technical Support

Please contact us with any technical questions about our products.

www.horiba.com/wq/support

User Support

- Our support website is available for registered customers and features:
- Data collection software
- Instruction manual downloadsMeasurement tips, etc.

www.horiba.co.jp/register

Validation Support

Please contact us with any questions or requirements for your validation procedure. •Traceability certification* •IQ/OQ/PQ support* •SOP guidance •FAQ *Optional services

Please read the operation manual before using this product to assure safe and proper handling of the product.

• The contents of this catalog are subject to change without prior notice, and without any subsequent liability to this company.

- The color of the actual products may differ from the color pictured in this catalog due to printing limitations.
- It is strictly forbidden to copy the content of this catalog in part or in full.
- All brand names, product names and service names in this catalog are trademarks or registered trademarks of their respective companies.
- Windows is a registered trademark of Microsoft Corporation in the United States and other countries

